CA[®] POWER GRIP Set

**Bedienungsanleitung
Operation Manual
Mode d'emploi
Manual de Instrucciones**

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Pablo Echarri, D.D.S.

The CA® POWER GRIPS were designed by Dr. Pablo Echarri, with the help of Julio Blasco drawing the first drafts and José Fernández making the first real prototypes, both of them working in Ladent, Barcelona. The next step in developing the products were made in SCHEU-DENTAL, specially by Christian Scheu and Peter Stückrad, that have worked intensively to make it happen.

The CA® POWER GRIPS were designed to improve movements with the CA® CLEAR ALIGNER technique in a more effective and rapid way, reducing the steps of treatment.

The main concept in the MFM is to apply traction forces (elastics) in extrusion and rotation treatments at home and during the nights, and combining these forces with CA® CLEAR ALIGNER & attachments treatment during the day.

The helix shape design of the MDI and DMI allow the CA® CLEAR ALIGNER to control the mesio-distal inclination of the teeth offering more retention of the appliance over the teeth achieving the movements more easily.

The FALA and FALI increase the pressure of the CA® CLEAR ALIGNER in the gingival area of the teeth making possible the torque correction.

The BUB and BUM increase the vertical dimension minimizing the occlusal interferences and improving the movements of the teeth in cross-bite position.

The forms BUAL and BUAS are very useful to increase the torque during the incisors intrusion with more effective forces.

They can also maintain the intrusion forces while the patient is not using the CA® CLEAR ALIGNER, permitting molar extrusion.

All the CA® POWER GRIPS requires a special protocol of lab set-up.

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CA® POWER GRIP Set

REF 5338.1

Content: 18 silicone forms (9 versions, 2 pcs. each); including support

CA® POWER GRIP forms enable fabrication of individual attachments for use in the CA® CLEAR ALIGNER technique.

Set contains nine different forms which can be filled with light curing composite to be directly applied on the prepared tooth.

Curing with LED lamp.

(We recommend using Light Bond - Sealant Resin [02-0033] and Light Bond Resin - Paste [02-0031] or Tetric EvoFlow®) www.smile-dental.de

CA® POWER GRIP, 2 silicone forms each

1 MFM Multifunctional Forced Movement

For extrusion, mesial and distal rotation for use with CA® splints and elastics.

REF 5330.1

2 MDI Mesial / Distal Inclination

For setting upright frontal teeth, mesial /distal inclination.

REF 5331.1

3 DMI Distal / Mesial Inclination

For setting upright frontal teeth; distal / mesial inclination.

REF 5332.1

4 FALA Force Application LABIAL

For protrusion, retrusion and torque control up to 3°. For fabrication of activation surfaces on the teeth in the gingival area with vestibular positioning - for teeth 11, 13, 21, 23, 33, 43.

REF 5333.1

5 FALI Force Application LINGUAL

For protrusion, retrusion and torque control up to 3°. For fabrication of activation surfaces on the teeth in the gingival area with labial, lingual and palatal positioning, for teeth 12, 22, 31, 32, 41, 42.

REF 5334.1

6 BUB Build-Up Bicuspid

Round pad for premolars (3 x 3 mm) for occlusal bite elevation (approx. 2 mm).

REF 5335.1

7 BUM Build-Up Molar

Oval pad for molars (3 x 5 mm) for occlusal bite elevation (approx. 2 mm).

REF 5336.1

8 BUAS Build-Up Anterior Small

Occlusal bite plane for intrusion of upper frontal teeth (2.5 x 2 mm)

REF 5337.1

9 BUAL Build-Up Anterior Large

Occlusal bite plane for intrusion of upper frontal teeth (2.5 x 4 mm)

REF 5339.1

CA® POWER GRIP Support made of stainless steel, sterilisable

REF 5340.1

2 pcs.

Product Description:

The CA® POWER GRIP technique enables fabrication of individual attachments.

Nine different forms are available which can be filled with light curing composite and directly applied on the tooth which has been prepared with Sealant/Resin. Then proceed with light curing.

Application:

1. Preparation of teeth: Polish the tooth surface by means of a rotating brush or with polishing paste. Then proceed with rinsing and air-drying the teeth.
2. Etching: Apply etching gel or liquid on a mixed paper. Dry the teeth to be etched and apply the etching agent on the corresponding tooth surfaces, avoiding rubbing and instead dabbing it. This takes 30 seconds, in case of milk teeth 120 seconds. Avoid any contact of soft tissue and etching gel.
3. Rinsing and drying: Rinse teeth thoroughly with water spray, each tooth for 10 seconds. In case of using etching gel for 20 seconds in order to stop etching and to remove dissolved enamel particles.
4. CA® POWER GRIP Form: Place the required CA® POWER GRIP form on the supplied support.
5. Application of primer: Apply light-curing sealant/liquid on the etched enamel surface, e.g. Assure or Light Bond Sealant.
6. Application of paste: Fill the clean form completely with light-curing primer or e.g. Light Bond past medium or LCR paste. If necessary, use a small instrument in order to fill undercuts completely. Positioning the attachment: Position the form on the required place for the attachment by pressing it slightly.
7. Curing of paste: Curing is done within approx. 5 - 20 seconds, depending on lamp intensity and primer.
8. Removal of form: After material has completely cured, proceed with carefully removing the form moving it slightly laterally.

Cleansing / Disinfection:

- The used form is to be cleansed by customary mild disinfection solutions.
- Make sure the form is completely clean before disinfection. Do not use any disinfection solutions containing formaldehydes or high-concentrated bleaching agents, as these may lead to damage and discolorations. Please follow the instructions of the corresponding manufacturer.

To be used only by expert staff!

The above information is for the doctor's use only, not for the patient. It is the doctor who has to evaluate this information with regard to each special treatment case and to inform the patient regarding advantages, potential side effects and complications.

HINTS:

This information corresponds to our actual experience to the best of our knowledge. However, this information is of non-binding nature and it is subject to any modification within the context of technical progress or development. This information only describes the nature and function of our products and our services, providing no warranty claim. The buyer is held responsible for a careful examination of the products' functions and applications by qualified expert staff. The buyer has to follow existing laws and regulations on his own responsibility.

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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©2013 Centro de Ortodoncia y ATM, Ladent, SL
All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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
1. MFM Multifunctional Forced Movements
2. MDI Mesial/Distal Inclination
3. DMI Distal/Mesial Inclination
4. FALA Force Application Labial
5. FALI Force Application Lingual
6. BUB Build-Up Bicuspid
7. BUM Build-Up Molar

The CA Power Grips were designed by Dr. Pablo Echarri, developed by Scheu-Dental, and tested by Ladent, SL and Echarri Clinic. The CA Power Grips consist of seven molds which allow easy attachment of composite retentions to the teeth, which facilitates certain dental movements with the CA Clear Aligner technique. The seven CA Power Grips are:

1. MFM Multifunctional Forced Movements
2. MDI Mesial/Distal Inclination
3. DMI Distal/Mesial Inclination
4. FALA Force Application Labial
5. FALI Force Application Lingual
6. BUB Build-Up Bicuspid
7. BUM Build-Up Molar

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In this Newsletter, we are going to talk about the CA Power Grip 4 (Force Application Labial - FALA) and CA Power Grip 5 (Force Application Lingual - FALI), used to correct the torque in approximately 5°.

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CA Power Grip 4 is bigger and it can be used in labial surface of upper incisors and upper and lower canines to increase the torque of these teeth.

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CA Power Grip 5 can be bonded in labial surface of upper lateral incisors and in labial surface of the four lower incisors, to increase the torque of these teeth. It can also be bonded in palatine surface from upper/lower canine to upper/lower canine, to reduce the torque of these teeth.

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Initial photographs of the case with the negative torque of lower central right incisor. The root of this tooth provokes the swelling of alveolar edge at this level.

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The impressions are taken and the torque of this tooth is increased in the set-up model. The torque is increased for 5°, controlled by CA Checker. Also, the lingual surface of this tooth is blocked with a Blokkout Resin layer, progressively thicker from incisal to gingival. If the impression is taken after the bonding of Power Grip, the model should be trimmed so CA Clear Aligner can carry out the pressure at this level.

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
This scheme shows how the Power Grip works, and how the pressure is carried out at the level of the Power Grip. The movement is possible because the space has been created when the CA Clear Aligner was made.

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*Isolate the working field with cotton rolls
separators and aspire.*

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*Etch with orthophosphoric acid during 20-
30 seconds.*

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**We recommend: Light Bond from Smile Dental, Tetric Eco Ceram from Ivoclar Vivadent, and Tetric Evo Flow from Ivoclar Vivadent.*

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*CA Power Grip 5 bonded on the tooth.
Frontal view.*

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*CA Power Grip 5 bonded on the tooth.
Lateral view.*

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*CA Power Grip 5 bonded on the tooth.
Occlusal view.*

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The alignment treatment in both arches using CA Clear Aligner has also been carried out. In these photographs, the final result can be seen. The upper midline matches with the center of a lower incisor, because the patient had only 3 incisors. The crown of upper right canine has been professionally replaced.

CLEAR-ALIGNER®

Antes/Before

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CLEAR-ALIGNER®

Antes/Before

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